1) DESCRIPTION AND DIRECTIONS FOR USE:

Avinent® Digital Health products are used to perform a patient-specific surgery/reconstruction during the planning phase and/or during reconstructive surgery.

The system consists of:

- **Anatomical models:** custom-made, single-use, acrylic or polyamide models designed to represent patient anatomy and used to improve and simplify surgery planning.
- **Surgical guides:** custom-made, single-use, acrylic or polyamide surgically-invasive guides designed and adapted to patient anatomy and used to improve, simplify and guide the surgical operation necessary to place titanium bone plates or other devices such as dental implants.
- **Bone plates:** custom-made, single-use, implantable, titanium alloy (Grade V - Ti6Al4V) plates designed and adapted to patient anatomy and planned surgical needs and used in the surgery to stabilize bone structures including bone grafts.
- **Customized prostheses:** custom-made, single-use, implantable, titanium alloy (Grade V - Ti6Al4V) or PEEK (polyetheretherketone) prostheses designed and adapted to patient anatomy and planned surgical needs and used in surgery to stabilize bone structures including bone grafts.
- **Bone screws:** series-made, single-use, implantable, titanium alloy (Grade V - Ti6Al4V), self-tapping, intraosseous fixation screws intended to attach bone plates and to stabilize and prevent movement of the bone segments / grafts.
- **Instruments:** series-made, reusable, stainless-steel twist-drills used to prepare bone to receive bone plates and screws during surgical procedures.
- **Accessories:** series-made, reusable screwdrivers, handles and instrument sterilization cassettes.

2) INSTRUCTIONS FOR USE

- Open the packaging
- The products listed below are supplied in a NON-STERILE condition. Clean and sterilise or disinfect the product before using it, as per point 5:
  - Anatomical Models
  - Surgical Guides
  - Bone plates
  - Customized prostheses
  - Bone Screws
  - Instruments
  - Accessories
- Were it to prove necessary to re-sterilise or to disinfect the product, we recommend you follow the process detailed in point 5 or that you contact the manufacturer.
- In the case of reusable products, clean and sterilise the product before using it, as per point 5:
- Use the product for the intended purpose.
- In the case of Surgical Guides, the steps you should follow to fit them on the patient are laid out below:
  - The guide is designed to adapt to the patient's anatomy.
  - The support surface must be completely free in order to ensure full adaptation.
  - The case report gives the position of the guide relative to the anatomy.
  - Do not apply too much pressure to the guide.
  - Ensure that you do not damage any critical anatomical structures during placement.
  - You may request anatomical models together with the guide, in order to test it before carrying out the surgery.
  - When you have the guide in a stable position, use the screws to fix it. Ensure that you have correctly identified the fixing holes.
  - If it is not possible to place the guide in the patient in a unique and stable position, the guide will not guarantee the precise transferral of the pre-op planning. Even in a stable position, it is possible that the guide cannot make contact with the bone along all its length, given that it is not always possible to determine the entire morphology. Undercuts depend on the shape of the patient's anatomy. When designing the guide, the undercut compensation is kept to a minimum in order to ensure maximum contact between the surface of the bone and the guide.
- During the cutting and the drilling it is recommended that you:
  - Ensure that you keep the guide in position on the contact surface while cutting and/or drilling.
  - Take all necessary precautions to avoid generating too much heat while cutting and/or drilling.
  - Do not use a saw blade that is thicker than the thickness indicated in the groove of the cutting guide.
  - Do not use a drilling tool that has a diameter that is larger than the one indicated in the hole. The case report specifies the drilling diameters that should be used.
  - Ensure that the saw blade follows the cutting groove or surface in order to obtain a correct osteotomy and to avoid cutting the groove.
  - Given that the inner diameters of the holes are larger than the diameter of the drill, try to drill along the central line in order to obtain a correct hole and to avoid perforating the inner wall of the hole.
- Use the screwdriver for the respective size of the bone screw.
- The surgeon is responsible for familiarising him or herself with the surgical techniques that are required for the implantation of Avinent® Digital Health products.

3) CONTRAINDICATIONS

- Infections of the surgical zone.
- Known allergies or hypersensitivity to the materials used when manufacturing the products.
- Patients who are rather unwilling to cooperate.
- Patients with serious uncontrolled systemic illnesses.

4) WARNINGS AND PRECAUTIONS

- Before using Avinent® Digital Health products, it is highly recommended that you should have been specifically trained in surgery techniques by a surgeon who is familiar with the Avinent® Digital Health.
- For optimum storage conditions, keep these products at room temperature and out of direct sunlight in a clean and dry place.
- Steam Sterilisation is not permitted in the case of models or guides manufactured in acrylic or thermolabile materials as to do so may give rise to prejudicial effects such as deformation and the loss of precision and integrity.
- It is the surgeon's responsibility to ensure that the products are completely sterile when they are used.
- Use specific products for cleaning, disinfection and sterilisation and validated procedures; regulate and inspect the devices used (disinfectant, steriliser) and ensure that the manufacturer's recommended parameters are respected during each cycle.
- AVINENT® Digital Health supplies all its articles with appropriate labelling that guarantees the correct traceability of the product in the patient's record, together with all the required documentation.
- Do not alter the products in any way. Altering the products may give rise to a setting that is inappropriate to the patient's anatomy. Were the products to be altered before or during surgery, the surgeon is solely responsible.
- Each product is correctly identified by a numerical code that corresponds to the batch and relates the product to the planned case.
- Before using the products, check that the numerical code or batch corresponds to the planned case.
- Reusing products intended to be used only once may create a risk of contamination and/or patient infection or may cause loss of mechanical integrity.
- Patients should be advised to consult the surgeon and seek medical attention if any movement of the plate or prosthesis is noted.
- Patients should be informed that special precaution should be taken before any exposure to MRT.
- In case of implant removal, the surgeon will have to proceed to unscrew the bone plate or customized prostheses.
- In case of exposed bone plate or customized prostheses, we recommend to hygiene them.
- Do not place dental implants in the mandible or maxilla until the healing process is finished. Patients should be informed that special precaution should be taken after dental implant placement.
- **Anatomic Models and Surgical Guides:**
  - Surgical models and guides are manufactured based on CAT images of the patient. If the anatomical characteristics of the patient have changed in any important way since the date of the CAT, the models should not be used. It is recommended that the model be used within 6 months from the date of the CAT.
  - Do not apply excessive force to the surgical models and guides and do not place heavy objects on them.
- **Instruments and Accessories:**
  - Do not store or transport instruments and/or accessories in such a way that they come into contact with each other, in order to avoid any possible damage to the product.
  - Do not use damaged instruments and/or accessories. It is recommended that you replace damaged instruments and/or accessories before they are used.
  - Clean the instruments and/or accessories as soon as possible after use, following the instructions given in point 5.

5) CLEANING AND STERILISATION INSTRUCTIONS

According to EN ISO 17664 standard, any medical device not marked as "sterile" must be cleaned, disinfected and/or sterilized before each use, in order to prevent any form of contamination. (This applies to the first use and, in the case of a reusable devices, all subsequent uses).

\* Instruments or equipments such as motors and contra-angle handpieces are excluded; in these cases it is recommend that you check each instrument's usage instructions.

GENERAL RECOMMENDATIONS	<ul style="list-style-type: none"> <li>• The user is responsible for the sterility of the product as well as for the damaged or dirty products.</li> <li>• Wear personal protective equipments (gloves, glasses, etc.)</li> <li>• Use disinfecting solutions which are approved for their efficacy.</li> <li>• The repetition of the indicated processes has a minimal effect on the products. The degradation of the products is due to the repetitive use.</li> </ul>
-------------------------	---

INSTRUCTIONS

Pre-disinfection:	<ul style="list-style-type: none"> <li>• In the case of reusable products, immerse all the products immediately after use in a disinfectant solution combined if possible with proteolytic enzymes. It is recommended that you ask your supplier for further details about disinfectant products.</li> <li>• Follow instructions of concentrations and immersion time given by the manufacturer (an excessive concentration may cause corrosion on the products).</li> <li>• Do not use any disinfectant solutions that contain any product that is not compatible with the products.</li> <li>• Rinse the products with purified water.</li> <li>• Should any impurities be visible on the products, a pre-cleaning is recommended by brushing the products manually.</li> </ul>
-------------------	---

Automatic cleaning and disinfection: (*)	<ul style="list-style-type: none"> <li>• Disassemble the products, if necessary.</li> <li>• Put them on an instrument rank and load it into the washer disinfectant.</li> <li>• Start the thermo disinfection cycle:                             <ul style="list-style-type: none"> <li>- Pre-cleaning with cold water 2 min</li> <li>- Cleaning 55°C (131°F) 5 min</li> <li>- Neutralizing with hot water 3 min</li> <li>- Intermediate rinsing with hot water 2 min</li> <li>- Final rinsing 93°C (199°F) 5 min</li> <li>- Drying 99°C (210°F) 35 min</li> </ul> </li> </ul> <p>The information provided is based on the use of:</p> <ul style="list-style-type: none"> <li>- Washer disinfectant: Miele G7883CD Program Vario TD</li> <li>- Chemicals: Neodisher FA (0.5%) as an alkaline detergent; ph 11.4-11.9 Neodisher N (0.1%) as a neutralizing agent; ph 2.21-1.4</li> </ul> <ul style="list-style-type: none"> <li>• Special instructions provided by the manufacturer of automated washing machines must be followed. Cleaning disinfectants are available commercially.</li> <li>• Discard any product with important defects.</li> <li>• Avoid any contact between products when placing in the thermo disinfectant.</li> <li>• Sufficient rinsing steps should be programmed.</li> <li>• Use "Purified Water" for rinsing.</li> <li>• Use filtered air for drying the products.</li> <li>• Use only approved thermo disinfectors according to EN ISO 15883.</li> </ul>
Manual cleaning and disinfection: (*)	<ul style="list-style-type: none"> <li>• Disassemble the products, if necessary.</li> <li>• Put them on an instrument rank.</li> <li>• Immerse them in a cleaning and disinfectant solution in an ultrasonic bath. Keep the products inside the ultrasonic bath for 15 min at 35°C (95°F) with 39K Hz applied. (To avoid protein coagulation Always keep a 1°&lt;45°C (113°F))</li> <li>• The solution has to be prepared following instructions given by the manufacturer.</li> </ul> <p>We recommend the following chemical:</p> <ul style="list-style-type: none"> <li>- Chemical: Instrunet Enzimático GPL (Inibsa) (0.83%)</li> </ul> <ul style="list-style-type: none"> <li>• Disinfection agents are available commercially. Please ask your supplier for details.</li> <li>• Rinse the products 2 times consecutively with purified water.</li> <li>• Dry the products with filtered air.</li> <li>• Discard any product with important defects.</li> <li>• Avoid any contact between products when placing in the solution.</li> </ul>
Cold disinfection (specifically for anatomic models and surgical guides manufactured in thermolabile acrylic material)	<ul style="list-style-type: none"> <li>• Place the products on a instrument rank or a perforated basket</li> <li>• Immerse the products in a bath filled with a disinfectant solution (e.g. Cidex Opa). The solution must be prepared in accordance with the manufacturer's instructions.</li> </ul>
Inspection:	<ul style="list-style-type: none"> <li>• Inspect the products and discard those with defects that can affect the resistance, safety or performance of them.</li> <li>• Assemble the devices, if necessary.</li> <li>• Dirty products must be cleaned and disinfected again.</li> </ul>
Packing:	<ul style="list-style-type: none"> <li>• Wrap the devices in sterilisation packaging.</li> <li>• Use FDA cleared commercially available wrap and packaging materials according to EN ISO 11607-1.</li> </ul>
Sterilisation: (*)	<ul style="list-style-type: none"> <li>• Sterilize the products by applying a fractionated pre-vacuum process according to EN 13060/ISO 17665-1.</li> <li>• Start the steam sterilisation cycle:                             <ul style="list-style-type: none"> <li>- 3 pre-vacuum phases with at 60 mbar minimum</li> <li>- Sterilization temperature: 134°C (273.2°F) / 135°C (274°F)</li> <li>- Minimum Holding Time: 4 minutes</li> <li>- Minimum Drying Time: minimum 14 minutes</li> <li>- Minimum Cooling Time: 10 minutes</li> </ul> </li> <li>• Respect maintenance of the autoclave device given by the manufacturer.</li> </ul>
Storage:	<ul style="list-style-type: none"> <li>• Keep the products in sterilisation packaging in a dry and clean environment.</li> <li>• Sterility cannot be guaranteed if packaging is open, damaged or wet.</li> </ul>

(\*) Not applicable to acrylics or thermolabile materials.

If you have any doubts, please consult our web: [www.avinent.com](http://www.avinent.com), send an email to [cmf@avinent.com](mailto:cmf@avinent.com) or call our customer helpline.

The AVINENT® Digital Health products have not been evaluated for safety and compatibility in the MR environment. The AVINENT® Digital Health products have not been tested for heating or migration in the MR environment.

6) SYMBOLS

As they appear on the product label:

 DO NOT REUSE LOT BATCH CODE MANUFACTURER <small>(For the EU only)</small> CE <small>CE MARKING without NB intervention - applicable for Class I, non-sterile products</small>	 NON-STERILE REF CATALOGUE NUMBER CONSULT INSTRUCTIONS FOR USE <small>(For the EU only)</small> CE <small>CE MARKING with intervention of NB 0197 - applicable for Class I/a / Class IIb products</small>
---	--

MANUFACTURER: AVINENT Implant System S.L. - Pol. Ind. Santa Anna I - E08251 Santpedor (Barcelona) - Spain

This document is subject to continuous revisions. Please verify that the current printed version is identical to the version that is available on the web: [www.avinent.com](http://www.avinent.com)

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